

LABORATORY EXPERIMENT

B.Sc. Semester: 7

PRACTICAL MANUAL

COURSE: (PMJDSC - Major)

(GROUP- A & GROUP- B)

COURSE CODE:SC23PMJDSCCHE701

GROUP –A

1. **INORGANIC QUALITATIVE ANALYSIS HAVING SIX REDICALS (MINIMUM EIGHT)**
 (OUT OF SIX ONE MUST BE FROM RARE GROUP AND ONE MUST BE
 TRACE IN A MIXTURE)

NEGATIVE REDICALS : Cl^- , Br^- , I^- , CO_3^{2-} , NO_3^- , NO_2^- , PO_4^{3-} , BO_3^{3-} , CrO_4^{2-}

POSSITIVE REDICALS : FROM 1ST AND 2ND GROUP AND VTH-B

RARE : FROM 1ST AND 2ND GROUP AND VTH-B

TRACE : Br^- , I^- , CO_3^{2-} , NO_3^- , NO_2^- .

Preliminary Test	Observation	Inference (Possible Ions)
Colour of Mixture	White	Na^+ , K^+ , NH_4^+ , Mg^{2+} , Ca^{2+} , Ba^{2+} , Al^{3+} , Zn^{2+} salts
	Blue	Cu^{2+}
	Yellow	Chromate (CrO_4^{2-}), Cd^{2+} salts
Solubility in Water	Completely soluble	May contain Na^+ , K^+ , NH_4^+ , NO_3^- , Cl^- , CH_3COO^- salts
	Partially soluble	Mixture may contain both soluble and insoluble salts
	Insoluble	May contain BaSO_4 , PbSO_4 , AgCl , CaCO_3 , SrSO_4 , metal oxides or sulphides

Dry Test

No.	Dry Test	Observation	Inference
1	Flame test	1. Green-blue flame	Cu^{+2} OR BO_3^{-3}
		2. Pale violate	K^{+1}
		3. Golden yellow	Na^{+1}
		4. Crimson red	Li^{+1}
2	Sodium hydroxide test Mix+NaOH	1. Ammoniacal Smell	NH_4^{+1}
		2. Yellow precipitate of mercury(II) oxide	Hg^{+1}
3	Mixture + dil. HCL (Heat mix. Slowly)	1. Break effervesces of CO_2	CO_3^{-2}
		2. After heating Yellow- red color gas evolved	NO_2^{-1}
		3. Add DW -White ppts	Pb^{+1} (1 st group may be Present)
4	Mixture + $\frac{1}{4}$ amount of MnO_2 + Con. H_2SO_4 (Heat mix. Slowly)	1. Violet color gas evolved	CO_3^{-2}
		2. Reddish-brown vapours are evolved which turn starch paper yellow	Br^-
		3. Reddish-brown vapours are evolved which turn starch paper yellow which don't turn starch paper yellow	NO_3^{-1}
		4. Greenish-yellow gas evolved gives white fumes with Ammoniacal solution	Cl^-
4	Mixture + Con. HCl (Heat mix. Slowly)+ Con. HNO_3 (excess amount)+ Ammonium Molibledom	1. Yellow ppt,s	PO_4^{-3}
5	In silica dish Mixture + Con. H_2SO_4 (Mix well) add some amount of Ethanol, Heat Mixture and fire it.	Green flame	BO_3^{-3}
6	Water extract (W.E.) + Zinc Urenileacetate	White Ppt,s	Na^+
7	Water extract (W.E.) + Picric acid	Yellow ppt's	K^+

8	Water extract (W.E.) [Yellow coloured + dil. HCL + Pass H ₂ S]	Yellow to Orange colour turns Green colour	CrO ₄ ⁻²
9	Soda extract + dil. H ₂ SO ₄ + small amount of Zn dust + heat mix add +KI (s) + Starch Solution	Blue Colour solution	NO ₃ ⁻¹ Presence (Br ⁻)

)

Wet Test for Positive Radical

Preparation of the Original solution:

Take small quantity of the Mixture and examine the solubility in the following solvents in the order given:

1. Water (Cold and hot)
2. Dil. HCl (Cold and hot)
3. Conc. HCl (Cold and hot)
4. Dil. HNO₃ (Cold and hot)
5. Conc. HNO₃ (Cold and hot)
6. Aqua Regia (3 parts of Conc. HCl + 1 part of Conc. HNO₃, mix and use the solution immediately, do not keep it at all)

Group Separation for Positive Radical

Original Solution (O.S.) + Add Dilute HCl

Precipitate → Group I

(Ag⁺, Pb²⁺, Hg⁺) & (Tl; w)



Filtrate + Add H₂S (g) in Dilute HCl (Acidic Medium) **Precipitate → Group II (Cu²⁺, Cd²⁺, Bi³⁺, As³⁺, Sb³⁺, Sn²⁺) & (Se;Te;Mo)**



Filtrate used for test of **(Na⁺, K⁺, NH₄⁺) and Li⁺**

Group Classification-1

The **Group I (Silver Group)** residue is washed using 2 cm³ of dilute HCl (see instruction 9) and then with 2 cm³ of water. Discard the washings. Boil residue in the centrifuge tube with 2 cm³ of water and centrifuge while hot. Decant the centrifugate into a test tube. Repeat extraction with 2 cm³ boiling water and combine the centrifugate.

Residue: White – Hg₂Cl₂ or AgCl. Yellow - H₂WO₄.
Boil the residue with 5 cm³ of dilute NH₃ solution and centrifuge.

Centrifugate: May contain Pb²⁺ or Tl⁺, which may crystallize on cooling. Carefully add 1 cm³ of conc. H₂SO₄ and evaporate till fumes are seen. Cool and dilute carefully to 5 cm³ with water. If a precipitate is formed, centrifuge.

Residue: Black precipitate, consisting of a mixture of Hg(NH₂)Cl and Hg.

Dissolve the residue by heating with 2 cm³ of aqua regia, cool and do the following tests:

1) To 1 cm³ of the solution, add NaOH solution in drops with shaking till excess. A black precipitate of Hg₂O is formed.

2) To 1 cm³ of the solution, add K₂CrO₄ solution and heat. A brown precipitate changing to red Hg₂CrO₄ on heating.

3) **Spot test:** Place a drop of faintly acid test solution (original) on a drop reaction paper and add a drop of saturated potassium nitrite (KNO₂) solution. A black spot of Hg is produced.

**Presence of Hg(I)
is confirmed.**

Centrifugate: May contain [Ag(NH₃)₂]⁺ or H₂WO₄. Add dilute HCl in drops with shaking till a faint precipitate forms, which is re-dissolved by adding one drop of dilute NH₃ solution. Add KI solution drop wise with shaking till precipitation, if any, is complete. Centrifuge.

Residue: Pale yellow precipitate of AgI.

1) To a portion of the above precipitate, add excess of sod. thiosulphate solution. The precipitate dissolves. Presence of Ag⁺ confirmed.

2) **Spot test:** Place a drop of the test solution (original) on a watch glass and add a drop of ammonium carbonate solution and mix. Withdraw a drop of clear liquid from the mixture and place it on a drop reaction paper and add a drop of K₂CrO₄ solution. A red ring is obtained due to formation of Ag₂CrO₄.

**Presence of Ag(I)
is confirmed.**

Centrifugate: Tungstic acid (H₂WO₄) solution.

Reduce to half its volume by evaporation, cool and do the following tests:

1) To 0.5 cm³ of the solution, add ferrous sulphate solution. A brown precipitate of ferrous tungstate is formed. Add 1 cm³ of dilute HCl and boil. The precipitate turns white and then yellow due to formation of tungstic acid.

2) **Spot test:** Place a drop of the test solution (original) on a drop reaction paper and add a drop of conc. HCl. A drop of saturated SnCl₂ solution is placed at the centre of the spot. A blue colour is produced due to formation of tungsten blue, W₂O₅.

**Presence of W(VI)
is confirmed.**

Residue: White – PbSO₄.
Dissolve in ammonium acetate solution (or a mixture of NH₃ solution and acetic acid) and the following tests are done:

1) To 1 cm³ of the solution, add acetic acid and K₂CrO₄ solution. A yellow precipitate of PbCrO₄ is formed.

2) To 1 cm³ of the solution, add KI solution. Boil and cool. A yellow precipitate of PbI₂ is formed, which dissolves on heating and re-precipitates as golden spangles on cooling.

3) **Spot test:** Place a drop of the solution on a drop reaction paper and add a drop of ammonium sulphide solution. A black spot of PbS is formed.

**Presence of Pb (II)
is confirmed.**

Centrifugate: May contain Tl⁺. Add dilute NH₃ solution and boil off excess NH₃. Cool and do the following tests:

1) To 1 cm³ of the solution, add acetic acid and K₂CrO₄ solution. A yellow precipitate of Tl₂CrO₄ is formed.

2) Flame test: Evaporate 2 cm³ of the solution to dryness in a china dish, cool and mix residue with a drop of conc. HCl using a glass rod. A speck of the mixture is taken at the tip of a nickel spatula and introduced into a Bunsen flame. An intense green colour appears

3) **Spot test:** Place a drop of the solution on a drop reaction paper and add a drop sodium thiosulphate solution. Add a drop of KI solution over it. A yellow spot of Tl (I) is produced.

**Presence of Tl (I)
is confirmed.**

Group Classification-2

Wash the **Group II** residue with 2 cm³ of a solution of ammonium chloride saturated with H₂S and centrifuge. Discard centrifugate. Transfer residue into a small beaker and add 5 cm³ of yellow ammonium sulphide (amm. polysulphide) solution. Heat to about 50 to 60°C and maintain at this temperature for 5 minutes with stirring. If there is a precipitate, cool and centrifuge.

Residue: Black – HgS, PbS, Bi₂S₃ or CuS. Yellow – CdS.

Analyse as given in Table 2A.

Centrifugate: Add conc. HCl dropwise till just acidic (do not add excess!) and warm while stirring. A fine white or pale yellow precipitate is only sulfur. If a flocculant precipitate is formed, centrifuge. If not, discard the solution.

Residue: Yellow – SnS₂, Se or As₂S₃. Orange red – Sb₂S₃.
Brown -Te, SnS or MoS₃.

Analyse as given in Table 2B.

Group Classification-2(A)

The **Group IIA (Copper Group)** residue is washed two times using 2 cm³ of distilled water and centrifuged. Discard washings. Transfer residue to a small beaker, add 10 cm³ of dilute nitric acid and boil for 5 minutes. If some solid material remains, cool and centrifuge.

Residue: Black – HgS. Dissolve by boiling with 2 cm³ of aqua regia. Do the following tests with this solution: 1) To 0.5 cm³ of the solution, add NaOH solution dropwise with shaking till excess. A brownish red precipitate changing to yellow HgO. 2) To 0.5 cm³ of the solution, add Na₂CO₃ solution dropwise and shake till effervescence just stops. Then add KI solution dropwise with shaking till excess. A red precipitate of HgI₂ is formed, which dissolves in excess reagent. 3) Spot test: Place a drop of the test solution on a spot plate. Add a drop of amm.thiocyanate solution and a drop of cobalt acetate solution. A blue colour is produced. Presence of Hg(II) is confirmed.	Centrifugate: May contain nitrate of Pb, Bi, Cu or Cd. Add dilute H₂SO₄ till precipitation (if any) is complete. Centrifuge.		
	Residue: White – PbSO₄. Dissolve in ammonium acetate solution (or a mixture of NH₃ solution and acetic acid) and the following tests are done: 1) To 1 cm³ of the solution, add acetic acid and K₂CrO₄ solution. A yellow precipitate of PbCrO₄ is formed. 2) To 1 cm³ of the solution, add KI solution. Boil and cool. A yellow precipitate of PbI₂ is formed, which dissolves on heating and re-precipitates as golden spangles on cooling. 3) Spot test: Place a drop of the solution on a drop reaction paper and add a drop of amm. Sulphide solution. A black spot of PbS is formed. Presence of Pb(II) is confirmed.	Centrifugate: May contain nitrate or sulphate of Bi, Cu or Cd. Add concentrated NH₃ solution dropwise with shaking till smell of ammonia persists. If a residue remains, centrifuge. Residue: White – Bi(OH)₃. The residue is divided into two parts. 1) To one part of the residue, add a few drops of conc. H₂O₂ solution. A yellowish brown precipitate of bismuthate ions (BiO₃⁻) is produced. 2) Dissolve the second part of the residue in 1 cm³ of dilute HNO₃ and divide into two parts. To one part, add KI solution dropwise with shaking to excess. A black precipitate of BiI₃ is first formed, which dissolves in excess KI to give an orange solution of [BiI₄]. 3) Spot test: Place a drop of the second part on a spot plate and add a drop of thiourea solution. An intense yellow colour is produced. Presence of Bi(III) is confirmed	Centrifugate: May contain [Cu(NH₃)₄]³⁺ or [Cd(NH₃)₄]³⁺. If deep blue in colour, Cu is present. Add dilute acetic acid with shaking till a clear solution is obtained. Divide into two parts. Divide the first part into three equal portions and do the following tests: 1) To one part, add pot. Ferrocyanide solution. A reddish brown precipitate of copper ferrocyanide is formed. 2) To the second portion, add pot.thiocyanate solution. A black precipitate turning slowly white is formed. $2 \text{ Cu(SCN)}_2 \rightarrow 2 \text{ CuSCN} + (\text{SCN})_2$. 3) Spot test: Place a drop of the third part on a spot plate and add a drop of tartaric acid solution, followed by a drop of NaOH solution. An intense blue colour is produced. Presence of Cu(II) is confirmed Divide the second part into three equal portions and do the following tests: 1) To one part, add NaOH solution. A white precipitate of Cd(OH)₂ is formed. 2) Through the second portion, pass H₂S. A yellow precipitate of CdS is obtained. 3) Spot test: Place a drop Cadion - 2 B reagent on a drop reaction paper. Place a drop from the third part over it, followed by a drop of KOH solution. A bright pink spot surrounded by a blue circle is produced. Presence of Cd(II) is confirmed.

Group Classification-2(B)

Transfer the Group IIB (Tin Group) residue into a small conical flask and place a funnel in its mouth. Add 5 cm ³ of conc. HCl and boil gently for 5 minutes. Dilute with 3 cm ³ of distilled water, cool and centrifuge.					
Residue: Yellow - Se or As ₂ S ₃ . Brown -Te or MoS ₃ . Dissolve in 5 cm ³ of conc. HCl containing a pinch of solid potassium chlorate (KClO ₃). Transfer solution to a small beaker and concentrate by evaporating on a water bath. Cool and add NH ₃ solution with stirring till smell of ammonia persists. Add 2 cm ³ of Mg(NO ₃) ₂ solution and stir for 5 minutes. When precipitation is complete, centrifuge.			Centrifugate: May contain Sb or Sn as chloride. Boil to expel H ₂ S. Cool and carry out the following tests using small portions of the solution.		
Residue: White crystalline Mg(NH ₄)AsO ₄ ·6H ₂ O Shake the residue well with 3cm ³ of water and divide into three portions. 1) To one portion, add AgNO ₃ solution containing a few drops of acetic acid. A red precipitate of silver arsenate, Ag ₃ AsO ₄ , is formed. 2) To the second portion, add 0.5 cm ³ of conc. HNO ₃ and excess of amm. molybdate solution. A yellow precipitate of amm. arsenomolybdate is formed. 3) Spot test: To the third portion, add two drops of conc. H ₂ O ₂ solution and warm. Acidify with a few drops of acetic acid. Place a drop of this mixture on a drop reaction paper followed by a drop of AgNO ₃ solution. A red spot is obtained. Presence of As is confirmed.	Centrifugate: May contain Se, Te or Mo as chlorides. Transfer to a small beaker and boil off NH ₃ . Add 2 cm ³ of conc. HCl and evaporate to half volume on a water bath. Add 2 cm ³ of a saturated solution of sodium sulphite (Na ₂ SO ₃) solution. If a precipitate is formed, centrifuge.		To 2 cm ³ of the solution in a test tube, add a pinch of iron filings and boil to reduce Sn(IV) to Sn(II). Cool and centrifuge. With draw the clear supernatant carry out the following tests. 1) to 1 cm ³ of the solution in a test tube, add two drops of mercuric chloride (HgCl ₂) solution. A silky white or grey precipitate of Hg ₂ Cl ₂ is obtained. 2) To 0.5 cm ³ of the solution in a test tube, add NaOH solution dropwise with shaking. A white precipitate of Sn(OH) ₂ soluble in excess NaOH is obtained. 3) Spot test: Place a drop of the solution on a drop reaction paper and add a drop of pot. Thiocyanate solution. Then add a drop of SnCl ₂ solution. A red spot of [Mo(SCN) ₆] ³⁻ is seen. Presence of Mo (VI) is confirmed.	To 1 cm ³ of the solution in a test tube, add just enough NH ₃ solution to neutralize. (A slight precipitation may occur). Add a pinch of solid oxalic acid and boil. Pass H ₂ S through the solution. An orange precipitate of Sb ₂ S ₃ is obtained. Spot test: Place a drop of the solution on a spot plate and add a minute crystal of sodium nitrite (NaNO ₂) . Then add a drop of rhodamine- B reagent. A blue colour is obtained. Presence of Sb is confirmed.	
	Residue: Red – Se. Dissolve the residue by heating with 1 cm ³ of conc. HNO ₃ . Neutralise by adding amm. carbonate solution dropwise till effervescence stops. Divide the solution into three parts. 1) To one part, add BaCl ₂ solution. A white precipitate of BaSeO ₃ is formed. 2) To the second part, add CuSO ₄ solution. A bluish green precipitate of CuSeO ₃ is formed. 3) Spot test: Place a drop of the third portion on a spot plate and add a drop of thiourea solution. A red precipitate of Se is formed. Presence of Se (IV) is confirmed.	Centrifugate: Dilute with equal volume of water. Add 1 cm ³ of KI solution and a drop of sodium sulphite (Na ₂ SO ₃) solution. If a precipitate is formed, centrifuge Residue: Black – Te Dissolve in 2 cm ³ of cold dilute HCl. Divide the solution into three parts. 1) Neutralise one portion by adding amm. carbonate solution dropwise till effervescence stops. Add BaCl ₂ solution. A white precipitate of BaTeO ₃ is formed. 2) To the second part, add excess of KI solution. A red colour is produced due to formation of [TeI ₆] ²⁻ . 3) Spot test: Place a drop of the third portion on a spot plate and add a drop of 50% hypophosphorus acid. Evaporate. A black stain due to Te is seen. Presence of Te (IV) is confirmed.			
	Presence of Sn is confirmed.				

Group Classification-5(B)

The Cations of this group are generally water soluble, hence they are detected from the water extract. Water is prepared by boiling 0.2g of the mixture with 1 cm³ of water. It is centrifuge after cooling and the centrifugate is used as follows.

For lithium

- 1) **Flame test:** To two drops of the test solution in a china dish, add a drop of conc. HCl and evaporate to dryness. Take a speck of the dry residue at the tip of a charred match stick and introduce into a colourless Bunsen flame near the base. The flame is coloured carmine red, indicating Li.
- 2) **Spot test:** Place a drop of the test solution on a spot plate. Add a drop of NaCl solution and a drop of ferric periodate reagent. A white precipitate of KLiFe[IO₆] is formed.

Presence of Li(I) is confirmed.

For Potassium

1. Take 5-6 drops of water extract and add 0.5 cm³ of 4N NaOH solution to it. Boil strongly to remove NH₄⁺ as NH₃ gas, cool the solution and used as follows

Acidify 2-3 drops of solution with 4N CH₃COOH and add 3-4 drops of 5% Na₂Co(NO₂)₆ solution to it. shake well and allow it to stand for about 10-15 sec – A yellow Crystalline Ppt

, K⁺ Present

C.T.

- a) 3 drops of the solution (free from NH₃) are treated with 4-5 drops of saturated picric acid solution. On shaking Well a - A yellow Crystalline Ppt
- b) 2 drops of the solution (free from NH₄⁺) + 3-4 drops of 50% EtOH + 4-5 drops of saturated tartaric acid solution – White Ppt

Spot Test : A drop of solution on a paper containing 1% Dipicryl amine reagent – orange red spot developed

For Sodium

Flame Test : Perform flame test with the mixture- Golden Yellow Flame

Na⁺ Confirm

SPECIAL TEST FOR SOME POSITIVE IONS

Rare Cation	Reagent / Spot Test	Observation	Inference
NH₄⁺ (Ammonium)	NaOH solution, warm	Pungent smell of NH ₃ ; moist red litmus turns blue	NH ₄ ⁺ present
Pb²⁺ (Lead)	Potassium chromate (K ₂ CrO ₄)	Yellow precipitate	Pb ²⁺ present
Cu²⁺ (Copper)	NH ₄ OH (excess)	Deep blue solution	Cu ²⁺ present
Bi³⁺ (Bismuth)	Thiourea in acidic medium	Yellow to orange colour	Bi ³⁺ present
Cd²⁺ (Cadmium)	H ₂ S in acidic medium	Bright yellow precipitate	Cd ²⁺ present
Hg²⁺ (Mercuric)	KI solution	Scarlet-red precipitate of HgI ₂ (soluble in excess KI)	Hg ²⁺ present
As³⁺ / As⁵⁺ (Arsenic)	H ₂ S in acidic medium	Yellow precipitate	Arsenic present
Sb³⁺ (Antimony)	H ₂ S in acidic medium	Orange precipitate	Antimony present
Sn²⁺ / Sn⁴⁺ (Tin)	HgCl ₂ solution	White precipitate turning grey/black	Tin present

SPECIAL TEST FOR SOME RARE POSITIVE IONS

Ion	Procedure (How to Perform the Test)	Observation	Inference
Tl⁺ (Thallium)	Take 2 mL original solution (O.S.) in a test tube. Add 2–3 drops of KI solution and shake.	Bright yellow precipitate of TlI forms, insoluble in excess KI.	Tl ⁺ present.
W (Tungsten, WO₄²⁻)	Take 2 mL O.S. Acidify with dilute HCl . Add 2–3 drops of freshly prepared SnCl₂ solution (or a small piece of zinc). Shake gently.	Deep blue colour (Tungsten Blue) appears.	Tungsten present.
Se (Selenium, SeO₃²⁻/SeO₄²⁻)	Take 2 mL O.S. Acidify with dilute HCl add fresh SnCl₂ solution and warm gently.	Brick-red precipitate of elemental selenium appears.	Selenium present.
Te (Tellurium, TeO₃²⁻/TeO₄²⁻)	Take 2 mL O.S. Acidify with dilute HCl and add fresh SnCl₂ solution . Warm gently if required.	Black precipitate of elemental tellurium appears.	Tellurium present.
Mo (Molybdenum, MoO₄²⁻)	Take 2 mL O.S. Add 1 mL dilute HCl , then 2–3 drops of NH₄SCN solution , followed by 2–3 drops of freshly prepared SnCl₂ solution . Mix well.	Deep blood-red colour develops due to thiocyanato-molybdenum complex.	Molybdenum present.
Li⁺ (Lithium)	Clean a platinum/nichrome wire by dipping it in concentrated HCl and heating until no colour is seen. Dip the wire into the sample and place it in the non-luminous Bunsen flame . Observe the flame colour.	Crimson (carmine) red flame .	Li ⁺ present.